





ANWASH
SURGICAL

DENTAL INSTRUMENTS



www.anwashsurgical.com

TABLE OF CONTENTS

Diagnostic	1
Probe	
Explorers	
Mirror Handles	
Mouth Mirrors	
Dressing Pliers	
Curettes	3
Gracey Standard	
Universal Curettes	
Periodontal Instruments (Scalers)	4
Restorative	5
Plastic Filling Carvers	
Plastic Filling (Titanium Nitride Coating)	
Chisels	
Excavators	
Condencers/Pluggers	
Burnishers	
Cravers	
Knives	
Waxing Instruments	
Margin Trimmers	
Endodontic Spreaders/Pluggers	
Restorative / Sinus Lift	9
Cord Packers	
Discoid-Cleoid	
Implant Instruments	10
Trephine Drills	
Osteotomes	
Surgery	11
Bone Splitting Instruments	
Periosteal Elevators	
Cheek/Lip Retractors	
Bur and File Stands	
Bone Management	12
Bone Files	
Curettes, Scrapers, Chisels, Gouges	
Restorative / Crown Removers	13
Matrix Retainers	
Syringes / Measurement	14
Endodontic / Restorative	15
Rubber Dams	
Wax, Plaster Knives, Cement Spatulas	
Impression Trays / Restorative	16
Amalgam Carriers Well & Gun	
Extracting Forceps	17
Forceps-English	
Forceps-Children	
Forceps-Special	
Forceps-Universal / American	
Pliers / Elevators	20
Orthodontic Pliers	
Root / Elevators	
Luxating Elevators	
Splinter Elevators	
Instruments Cassettes	22

ANWASH
SURGICAL




DIAGNOSTIC




PROBES

HANDLE OPTIONS



Available with/without
colour marking

EXPLORERS



Probes are precision fabricated to provide accuracy for intraoral gauging. Probes with extra fine, sharp points for diagnosing tartar or concretions, caries, filling and tooth morphology. All probe instruments are designated in millimeters. Avoid temperature above 200°C / 398°F, at which point your Hi-Lite Probe markings may become discolored.

HANDLES FOR MIRRORS

Handles with M2.5 Thread



MOUTH MIRRORS

Note: All Mouth Mirrors are Available in Plane & Magnifying Glass




Pack of 12

AS-MM126
Back Surface


AS-MM127
Rhodium Front Surface



Front surface is rhodium coated, reflective coating is on the upper side of the mirror glass, in this way a clear and undistorted image is achieved.



DRESSING PLIERS



CURETTES

GRACEY STANDARD

HANDLE OPTIONS



Silicone Rubber 9.5 mm



Colour Resin 9.5 mm



Stainless Steel 9.5 mm



AS-G134
Gracey 1-2



AS-G135
Gracey 7-8



AS-G136
Gracey 13-14



AS-G137
Gracey 3-4



AS-G138
Gracey 9-10



AS-G139
Gracey 15-16



AS-G140
Gracey 5-6



AS-G141
Gracey 11-12



AS-G142
Gracey 17-18



Gracey 5-6

Application field of original
Gracey Curettes
1-2, 3-4, 5-6
Front teeth

7-8, 9-10
Bicuspid and molars

11-12, 15-16
Bicuspid and molars mesial,
furcation: mesial
part of the distal root

13-14, 17-18
Bicuspid and molars distal,
furcation: distal part of the mesial
root.

UNIVERSAL CURETTES



8 mm



8 mm

HANDLE OPTIONS



AS-UC143
Barnhart 1-2



AS-UC144
McCall 13-14



AS-UC145
McCall 17-18 S



AS-UC146
Columbia 2L-2R



AS-UC147
Langer 1-2



AS-UC148
Barnhart 5-6



AS-UC149
McCall 13-14 S



AS-UC150
Columbia 4L-4R



AS-UC151
Langer 3-4




AS-UC152
Rules 3-4



AS-UC153
McCall 17-18



AS-UC154
Younger-Good 7-8



AS-UC155
Langer 5-6



Columbia 4L-4R



SCALERS

HANDLE OPTIONS



6. mm



6 mm



9.5 mm



6.5 mm



AS-S156
Jacquette 3-4



AS-S159
Jacquette 14-15



AS-S162
Jacquette 3-U5



AS-S157
Jacquette 3-U15



AS-S160
Jacquette 4-U5



AS-S163
Jacquette 4-U15



AS-S158
M23



AS-S161
U15-30



AS-S164
U15-33

8 mm

H 6 -7



AS-S165
Eagle Claw



AS-S166
Curaler - N 137



AS-S167
S 311-312



AS-S168
N3



AS-S169
Morse 0-00



AS-S170
Eagle Talon



AS-S171
N135



AS-S172
N1



AS-S173
SE Nebraska 128
(Single End)



AS-S174
N2


SurgiDoc instrument is crafted by hand. Tips are fabricated individually by highly skilled precision technicians. Each instrument is continually held to strict quality controls during each stage of manufacturing, with tapering point working end. Well suitable for inter-dental cavities.

PLASTIC FILLING CARVERS



PLASTIC FILLING

TITANIUM NITRIDE COATING IS SIGNIFICANTLY HARDER THAN STEEL.



CHISELS



6 mm Stainless Steel Double End (Standard)



SD-C191
Wedelsteadt 3-4
(11.5-15-3)



SD-C192
Binangle 8-9
(20-9-8)




SD-C193
Wedelsteadt 5-6
(15-15-3)



SD-C194
Binangle 11-12
(15-8-8)

SurgiDoc excavators, used for the removal of softened carious dentin, are equipped with spoon-shaped working end for work on deep areas and a round working end for tissue-saving cavity preparation.

EXCAVATORS




6 mm




SD-EX195
17



SD-EX196
20




SD-EX197
42-43
(15-8-14)




SD-EX198
70-71




SD-EX199
E-1
(12-9-15)



SD-EX200
18




SD-EX201
38-39
(11.5-7-14)




6.6 mm




SD-EX202
64-65




SD-EX203
74-75



SD-EX204
E-2
(15-9-15)



SD-EX205
40-41
(11-7-14)




SD-EX206
66-67



SD-EX207
78-79



SD-EX208
E-3
(25-9-15)




8 mm


74-75

CONDENSERS / PLUGGERS (NON-SERRATED)


HANDLE OPTIONS



BURNISHERS



CARVERS



KNIVES



WAXING INSTRUMENTS




MARGIN TRIMMERS



ENDODONTICS SPREADERS / PLUGGERS

Available with/without colour marking




RESTORATIVE / SINUS LIFT

CORD PACKERS





DISCOID-CLEOID




All of our instruments are hand sharpened by highly trained technicians to provide sharper, longer lasting tips and edges.

SINUS LIFT ELEVATORS




TREPHINE DRILLS




OSTEOTOMES

Handle Size: 10 mm, 13 mm



BONE SPLITTING INSTRUMENTS




PERIOSTEALS ELEVATOR

For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue.



SD-P300
Elevator / Periodontal Chisel



SD-P303
Glickmann



SD-P301
Buser / Dissector



SD-P304
Elevator Papilas



SD-P302
Buser / Papilla Remover



SD-P305
PPLABANCA 6



SD-P307
Williger
Periosteal Elevator
4,5 mm. 6.0 mm

CHEEK, LIP RETRACTOR AND STANDS



SD-R308



SD-R309



SD-H310
Bur Holder



SD-H311
Endodontic Holder for Files and Reamers

BONE FILES




Fig. 1




Fig. 2




Fig. 3



AS-BF312
Bone Files
Miller- Colburn 1 X- Cut
18 cm




Fig. 4




Fig. 5

CURETTES, SCRAPERS, CHISELS AND GOUGES



AS-321
Miller & Lucas
Spoon Shape
2.2mm, 2.8 mm
6

AS-322
Will
2.5mm



AS-313
Articulating Paper Forcep
Miller
15 cm



AS-314
Williger 14.5 cm

Sharp Spoon



AS-315
Hemingwa
16.5 cm



AS-316
Will Sharp Spoon 16.5 cm

2,5 mm

3,5 mm



AS-317
Periosteal Elevator Freer (Sharp / Blunt)



AS-318
Bone Chisel Partsch 17 cm (3-6 mm)



Partsch **AS-319**
Gouge Width 4 mm, 5 mm, 6 mm, 7 mm

17 cm (3-6 mm)



AS-320
Gouge
17 cm (3-6 mm)

RESTORATIVE / CROWN REMOVERS

MATRIX RETAINERS



SD-MR323
Ivory
65 mm



SD-MR324
Tofflemire-Universal
60 mm



SD-MR325
Tofflemire-Junior
45 mm



SD-MR326
Tofflemire-Senior
50 mm



SD-MR327
Nystrom, Right Left
75 mm

CROWN INSTRUMENTS



SD-CI328
Orthodontic Band Pusher



SD-CI329
Crown Remover
Automatic with 5 insets including lock and key



SD-CI330
Crown-Butler with 3 insets
3/4, chrome plated brass



SD-CI331
Morrel Crown Remover
(Manual) with 5 insets

WATER AND CARTRIDGE SYRINGES



AS-S332

Aspirating Cartridge Syringe
metric / with US thread
screw tops for cannulas ointment cannula
blunt/curved
1.8 cc, 2.2 cc



AS-S333

Water Syringe
10 ml with Luer
lock cannula



AS-S334

Cartridge Syringe
(TKC), cylinder
without aspiration
metric / US threads
chrome plated brass
1.8 cc, 2.2 cc

Also Available With Aspiration



AS-S335

Skypject Intra lingual Syringe
Used in performing PDL
injections. Brass
and stainless steel
11 cm



AS-S336

Self Aspirating Syringe
anodized aluminum,
and stainless steel
1.8 cc, 2.2 cc

MEASURING INSTRUMENTS



AS-MI337

Positive Tube Holder
Anodized
Aluminium



AS-MI338

Boley Caliper
with Lock 14 cm
Scale 0-100 mm



AS-MI339

Measurement of
bone thickness with a precision
of up to 0.05 mm
14.5 cm



AS-MI340

Ivanson
Wax/Metal 10 cm (Scale 0-10)
New Dental Caliper





AS-MI341

Castroviejo Caliper
10 cm (Scale 0-40)

RUBBER DAM INSTRUMENTS

RUBBER DAM CLAMPS AS-C345



AS-F345-A
Rubber dam frame
big/small 15X15 cm
Stainless Steel/Plastic



AS-RD342
Rubber dam punch
Anisworth 16 cm
also available Ivory with 1 and 2 rows of holes



AS-RD343
Rubber dam
clamp forcep Brewer
17 cm



AS-RD344
Rubber dam clamp forcep
Stokes
17 cm

WAX AND PLASTER KNIVES, CEMENT SPATULAS

round 8.0 mm round 6.5mm round 5.0mm angular 8.0mm angular 6.5 mm angular 5.0 mm



AS-K346
Wax knives with
wooden handles
Fahenstock
13 cm,
17.5 cm



AS-K347
Lessman
with grooved blade
13 cm, 17 cm



AS-K348
Fahenstock
All metal wax knife
13 cm



AS-CP349



AS-CP350



AS-CP351



AS-CP352



AS-CP353



AS-CP354



AS-CP355

IMPRESSION TRAYS

Assorted Colours Moulded Plastic



AS-IT356

AS-IT357

AS-IT358

		mm	mm
U0	xxs	63	48
U1	xs	65	52
U2	s	68	54
U3	M	71	58
U4	L	73	59
U5	XL	78	63

Perforated Without Retention Rim



AS-IT359

AS-IT360

Perforated With Retention Rim



AS-IT361

AS-IT362

		mm	mm
L0	xxs	56	44
L1	xs	66	50
L2	s	70	52
L3	M	73	56
L4	L	76	58
L5	XL	78	60

Un perforated With Retention Rim



AS-IT363

AS-IT364

Also Available Partial Trays



AS-IT365
Perforated for
edentulous upper jaws

AS-IT366
Un perforated For
toothed lower jaws

		mm	mm
L		79	70
M		76	67
S		73	60

		mm	mm
L		79	70
M		76	67
S		73	60



AS-IT367
Perforated
For Partial
Impressions right
R1, R2, R3



AS-IT368
Unperforated
For Partial
Impressions Left
L1, L2, L3

AMALGUM CARRIER, WELL AND GUN



AS-AC369
23 cm

AS-AC370
Ivory 23 cm
Also Available Single Ended 17.5 cm

AS-AW371
Amalgam Well
H 1.5 cm / D 3.5 cm
H 3.5 cm / D 3.5 cm



13 cm, 15 cm
AS-AG372
Amalgum Gun
30° 45° 90°

ENGLISH PATTERN



AS-E373
Fig. 1
For Upper Incisors
and Canines, Wide



AS-E374
Fig. 168
For Upper Canines



AS-E375
Fig. 136
For Upper Premolars



AS-E376
Fig. 17
For Upper Molars, Right



AS-E377
Fig. 18
For Upper Molars, Left



AS-E378
Fig. 7
For Upper Premolars



AS-E379
Fig. 67
For Upper Wisdom Teeth, Narrow
Available: For Upper Wisdom Teeth, Wide



AS-E380
Fig. 29
For Upper Front Roots



AS-E381
Fig. 30
For Upper Roots



FORCEPS FOR CHILDREN



SD-E385
Fig. 89
For Upper Molars,
Right, Gripping In Depth



SD-E382
Fig. 18 A
For Upper Molars, Either Side



SD-E383
Fig. 51 A
For Upper Front Roots,
Narrow Beaks



SD-E384
Fig. 13
For Lower Canines And Premolars



SD-E386
Fig. 33 A
For Lower Roots, With
Narrow Beaks



SD-E387
Fig. 33
For Lower Incisors,
Premolars And Roots



SD-E388
Fig. 75
For Lower Canines
And Premolars



SD-E389
Fig. 24
For Lower Molars, Left



SD-E390
Fig. 31
For Lower Premolars
And Roots



SD-E391
Fig. 73
For Lower Molars



SD-E392
Fig. 21
For Lower Wisdom
Teeth And Molars

EXTRACTION FORCEPS

ENGLISH PATTERN FOR CHILDREN



AS-E393
Fig. 37
For Upper Incisors



AS-E394
Fig. 158
For Upper Molars,
Either Side, Small



AS-E395
Klein
For Upper Incisors



AS-E396
Schmidt
For Upper Roots



**Also Available Upper
And Lower Molars**



AS-E397
Klein
For Lower Incisors



AS-E398
Klein
For Lower Roots



AS-E399
Cryer Fig. 150 S
Tooth Forceps For Children
For Upper Incisors, Canines
And Premolars



AS-E400
Cryer Fig. 151 S
Tooth Forceps For Children
For Lower Incisors, Canines
And Premolars And Molars



AS-E406
Cryer Fig. 151
For Lower Incisors, Canines,
Premolars And Roots

AMERICAN AND UNIVERSAL PATTERN



AS-E401
Harris Fig. 18 R
For Upper Molars, Right



AS-E402
Fig. 203
For Lower Incisors,
Premolars And Roots



AS-E403
Fig. 17
For Lower Molars



AS-E404
Fig. 222
For Lower Wisdom Teeth



AS-E405
Hull Fig. 101
Universal Tooth
Forceps

ORTHODONTIC PLIERS



AS-OP407
Angle (Similar Instruments: Adams)
12.5 cm
Wire Bending
Plier



AS-OP408
Aderer max
12 cm, 13 cm



AS-OP409
Wire Cutter 12.5 cm
Tip : Tungsten Carbide
It cuts ligature wires and pins of any size
Cutting Capacity: archwires up to 0,56x0,71 mm



AS-OP410
Distal End Cutter 12.5 cm-safety hold
designed for distal arch cutting
in the patient's mouth.

SurgiDoc welcomes the opportunity to manufacture any other orthodontic pliers. Please send us actual sample and/or picture with detailed written description that includes name and measurements.

ROOT ELEVATORS



AS-RE411
Apical



AS-RE412
Apical
Fig. 27



AS-RE413
Apical
Fig. 28



AS-RE414
Apical
Fig. 301



AS-RE415
Apical
Fig. 302



AS-RE416
Apical
Fig. 303



AS-RE417
Apical
Fig. 304



AS-RE418
Kopp



AS-RE419
Kopp
Fig.3



D-RE4220
Kopp
Fig.4



AS-RE421
Kopp
Fig.3




AS-RE422
Kopp
Fig.4

ELEVATORS

ROOT AND LUXATING ELEVATORS

Luxating Instruments
Luxating instruments have thin, sharp blades that are used to cut periodontal ligaments.



SPLINTER ELEVATOR



ROOT ELEVATORS



Available in Blue, Red, Yellow and Green Silicon



Removable Rack
Twist Lock System



Removable/Fixed Rack
Twist Lock System



ANWASH SURGICAL



Roras Road, Sialkot - Pakistan



+92 310 4719254



info@anwashsurgical.com



www.anwashsurgical.com

